


# Leuk●Stain

For quantitative detection of leukocyte concentration in Mammalian semen samples



## Features:

- Rapid reaction time (<10 seconds).
- Minimal bubble formation.
- Very simple to use.
- Less operational time than our competition.

Detailed Product Description	
Method:	Peroxidase Staining
Storage Temperature:	2 - 8°C
Package:	1500 Tests/kit
Shelf Life:	24 Months

## General Sources regarding technology

- 1. World Health Organization. 1999. WHO laboratory manual for the examination of human semen and sperm-cervical mucus interaction. Cambridge University Press. Fourth Edition.
- 2. Tremellen K. Oxidative stress and male infertility--a clinical perspective. 2008. Hum Reprod Update. 4:243-258

## For more information

Email: [info@pylon-healthcare.com](mailto:info@pylon-healthcare.com) Website: [www.Pylon-Healthcare.com](http://www.Pylon-Healthcare.com) Phone: +20 122 555 6925, +20 1000 922 377 , +20 38520295

