



INFINITY MEDIA
YOUR DAILY CULTURE NEEDS

HEPES Buffer Solution (1X)

A BUFFERING SYSTEM TO MAINTAIN RELEVANT PHYSIOLOGICAL PH RANGE
AND OSMOLAR BALANCE IN CELL CULTURE APPLICATIONS.

LIQUID, STERILE FILTERED, IVF STANDARD FORMULATION
WITH PHENOL RED

SUITABLE FOR

- ➔ FOR OUTSIDE INCUBATORS DAILY GAMETE HANDLING
- ➔ GAMETES HANDLING, WASHING AND MANIPULATION

Osmolality	262 - 315 mOs/kg
pH	7.2 - 7.6
Antibiotic concentration	10ul/ml
Sterility	Pass
Salt Toxicity Test	Pass
Packaging	500ml sterile PET airsealed bottles

Shelf life recommended: 9months from production or 6 weeks from oppening

Website WWW.Pylon-Healthcare.com

Email info@Pylon-HEalthcare.com

PYLON
HEALTHCARE