



INFINITY MEDIA
YOUR DAILY CULTURE NEEDS

Earle's Balanced Salts Solution

A BUFFERING SYSTEM TO MAINTAIN RELEVANT PHYSIOLOGICAL PH RANGE AND OSMOLAR BALANCE IN CELL CULTURE APPLICATIONS.

LIQUID, STERILE FILTERED, SUITABLE FOR CELL CULTURE WITH SODIUM BICARBONATE, MAGNESIUM, CALCIUM & PHENOL RED

SUITABLE FOR

- ➔ AS A DILUENT, TRANSFER MEDIUM OR FOR SAMPLE WASHING.
- ➔ OOCYTE COLLECTION & INCUBATION POST COLLECTION
- ➔ ALL EUKARYOTIC CELL LINES DAILY CELL & TISSUE CULTURE

Osmolality	272 - 315 mOs/kg
pH	7.2 - 7.6
Glucose Concentration	0.9 - 1.1 g/l
Antibiotic concentration	10ul/ml
Sterility	Pass
Salt Toxicity Test	Pass
Packaging	500ml sterile PET airsealed bottles

Shelf life 1 year from production or 6 weeks from oppening

Website WWW.Pylon-Healthcare.com

Email info@Pylon-HEalthcare.com

PYLON
HEALTHCARE